

Select BioProducts Horizontal Electrophoresis System Accessories

A wide assortment of gel trays and combs are available for the Select BioProducts Horizontal Electrophoresis Systems. All the gel trays are UV transparent and work with the casting dams supplied with the original unit. Additional casting dams are available. Slots on the sides of the gel trays hold the combs in place. The combs consist of two pieces, allowing the depth of the comb within the gel to be adjusted. Many of the combs are multichannel compatible.

Mini 7 and Mini 10 Systems (7 cm wide)

Catalog No.	Description	Well vol.*
SBE1007-7-GT	Gel tray, UV transparent, 7 x 7 cm	NA
SBE1007-10-GT	Gel tray, UV transparent, 7 x 10 cm	NA
SBE1007-CD	Casting dams for 7 cm wide gel trays, pk of 2	NA
SBE1107-16-1.5	Comb, 16 teeth, 1.5 mm thick	15 µl
SBE1107-12MC1.5	Comb, 12 teeth, 1.5 mm thick, MC compatible	20 µl
SBE1107-10-1.5	Comb, 10 teeth, 1.5 mm thick	27 µl
SBE1107-5-1.5	Comb, 5 teeth, 1.5 mm thick	68 µl
SBE1107-16-1	Comb, 16 teeth, 1 mm thick	10 µl
SBE1107-12MC-1	Comb, 12 teeth, 1 mm thick, MC compatible	14 µl
SBE1107-8-1	Comb 8 teeth, 1 mm thick	25 µl

Midi 10 Systems (10 cm wide)

Catalog No.	Description	Well vol.*
SBE1010-10-GT	Gel tray, UV transparent, 10 x 10 cm	NA
SBE1010-CD	Casting dams, for 10 cm wide gel tray, pk of 2	NA
SBE1110-25-2	Comb, 25 teeth, 2 mm thick	20 µl
SBE1110-16-2	Comb, 16 teeth, 2 mm thick	32 µl
SBE1110-20MC1.5	Comb, 20 teeth, 1.5 mm thick, MC compatible	20 µl
SBE1110-16-1.5	Comb, 16 teeth, 1.5 mm thick	24 µl
SBE1110-12-1.5	Comb, 12 teeth, 1.5 mm thick	34 µl
SBE1110-10MC1.5	Comb, 10 teeth, 1.5 mm thick, MC compatible	41 µl
SBE1110-20MC-1	Comb, 20 teeth, 1 mm thick, MC compatible	14 µl
SBE1110-10MC-1	Comb, 10 teeth, 1 mm thick, MC compatible	27 µl
SBE1110-12-1	Comb, 12 teeth, 1 mm thick	23 µl
SBE1110-8-1	Comb, 8 teeth, 1 mm thick	41 µl
SBE1110-16-0.75	Comb, 16 teeth, 0.75 mm thick	12 µl
SBE1110-12-0.75	Comb, 12 teeth, 0.75 mm thick	17 µl


Midi Plus 10 and Midi Plus 15 Systems (15 cm wide)

Catalog No.	Description	Well vol.*
SBE1015-10-GT	Gel tray, UV transparent, 15 x 10 cm	NA
SBE1015-15-GT	Gel tray, UV transparent, 15 x 15 cm	NA
SBE1015-CD	Casting dams, for 15 cm wide gel tray, pk 2	NA
SBE1115-20-2	Comb, 20 teeth, 2 mm thick	43 µl
SBE1115-16MC-2	Comb, 16 teeth, 2 mm thick, MC compatible	54 µl
SBE1115-14MC-2	Comb, 14 teeth, 2 mm thick, MC compatible	49 µl
SBE1115-35-1.5	Comb, 35 teeth, 1.5 mm thick	15 µl
SBE1115-30MC1.5	Comb, 30 teeth, 1.5 mm thick, MC compat.	19 µl
SBE1115-20-1.5	Comb, 20 teeth, 1.5 mm thick	32 µl
SBE1115-18MC1.5	Comb, 18 teeth, 1.5 mm thick, MC compat.	17 µl
SBE1115-16MC1.5	Comb, 16 teeth, 1.5 mm thick, MC compat.	41 µl
SBE1115-14MC1.5	Comb, 14 teeth, 1.5 mm thick, MC compat.	44 µl
SBE1115-12-1.5	Comb, 12 teeth, 1.5 mm thick	61 µl
SBE1115-10-1.5	Comb, 10 teeth, 1.5 mm thick	68 µl
SBE1115-8-1.5	Comb, 8 teeth, 1.5 mm thick	100 µl
SBE1115-30MC-1	Comb, 30 teeth, 1 mm thick, MC compatible	13 µl
SBE1115-28MC-1	Comb, 28 teeth, 1 mm thick, MC compatible	11 µl
SBE1115-20-1	Comb, 20 teeth, 1 mm thick	21 µl
SBE1115-18MC-1	Comb, 18 teeth, 1 mm thick, MC compatible	11 µl
SBE1115-16MC-1	Comb, 16 teeth, 1 mm thick, MC compatible	27 µl
SBE1115-16-1	Comb, 16 teeth, 1 mm thick	27 µl
SBE1115-12-1	Comb, 12 teeth, 1 mm thick	44 µl
SBE1115-10-1	Comb, 10 teeth, 1 mm thick	29 µl
SBE1115-28MC.75	Comb, 28 teeth, 0.75 mm thick, MC compat.	8 µl

Maxi 20 System (20 cm wide)

Catalog No.	Description	Well vol.*
SBE1020-20-GT	Gel tray, UV transparent, 20 x 20 cm	NA
SBE1020-CD	Casting dams, for 20 cm wide gel tray, pk/2	NA
SBE1120-50-2	Comb, 50 teeth, 2 mm thick	21 µl
SBE1120-40MC-2	Comb, 40 teeth, 2 mm thick, MC compatible	23 µl
SBE1120-40MC1.5	Comb, 40 teeth, 1.5 mm thick, MC compatible	17 µl
SBE1120-20MC1.5	Comb, 20 teeth, 1.5 mm thick, MC compatible	41 µl
SBE1120-40MC-1	Comb, 40 teeth, 1 mm thick, MC compatible	11 µl
SBE1120-30-1	Comb, 30 teeth, 1 mm thick	17 µl
SBE1120-25-1	Comb, 25 teeth, 1 mm thick	21 µl
SBE1120-20MC-1	Comb, 20 teeth, 1 mm thick, MC compatible	27 µl
SBE1120-20MC.75	Comb, 20 teeth, 0.75 mm thick, MC compat.	20 µl
SBE1120-16-0.75	Comb, 16 teeth, 0.75 mm thick	30 µl

*Approximate sample volume for a 5 mm thick gel



10 x 10 cm casting tray with
10 tooth combs and end dams

Select BioProducts Mini PAGE Electrophoresis System Accessories

A variety of accessories are available to customize the Select BioProducts Mini PAGE Electrophoresis System to the needs of your laboratory. All modules fit in the running tank supplied with the original system. Combs, glass plates and spacers can be used with the original PAGE system.

PAGE Accessories (for 10 x 10 cm system)

Catalog No.	Description
SBE2010-PM	PAGE module (without plates, combs, caster)
SBE2110-NG-2	Notched glass plates, 2 mm thick, 2 ea
SBE2110-PG-2	Plain glass plates, 2 mm thick, 2 ea
SBE2110PG0.75BS	Plain glass plates w/.75 mm bonded spacers, 2 ea
SBE2110-PG-1-BS	Plain glass plates w/1 mm bonded spacers, 2 ea
SBE2110-PG1.5BS	Plain glass plates w/1.5 mm bonded spacers, 2 ea
SBE2110-PG-2-BS	Plain glass plates w/2 mm bonded spacers, 2 ea
SBE2110-NG.75BS	Notched glass plates w/.75 mm bonded spacers, 2 ea
SBE2110-NG-1-BS	Notched glass plates w/1 mm bonded spacers, 2 ea
SBE2110-NG1.5BS	Notched glass plates w/1.5 mm bonded spacers, 2 ea
SBE2110-NG-2-BS	Notched glass plates w/2 mm bonded spacers, 2 ea
SBE2110-DP	Dummy plate
SBE2110-0.75-S	Spacers, 0.75 mm thick, pack of 2
SBE2110-1-S	Spacers, 1 mm thick, pack of 2
SBE2110-1.5-S	Spacers, 1.5 mm thick, pack of 2
SBE2110-2-S	Spacers, 2 mm thick, pack of 2
SBE2110-12-1	Comb, 12 sample, 1 mm thick
SBE2110-12-1.5	Comb, 12 sample, 1.5 mm thick
SBE2110-10-1	Comb, 10 sample, 1 mm thick
SBE2110-10-1.5	Comb, 10 sample, 1.5 mm thick
SBE2110-10-0.75	Comb, 10 sample, 0.75 mm thick
SBE2110-5-1.5	Comb, 5 sample, 1.5 mm thick
SBE2110-5-0.75	Comb, 5 sample, 0.75 mm thick
SBE2110-16MC-1	Comb, 16 sample, 1 mm thick, MC compatible
SBE2110-20-1	Comb, 20 sample, 1 mm thick
SBE2110-5-1	Comb, 5 sample, 1 mm thick
SBE2110-20-0.75	Comb, 20 sample, 0.75 mm thick
SBE2110-12-2	Comb, 12 sample, 2 mm thick
SBE2110-8MC-1	Comb, 8 sample, 1 mm thick, MC compatible
SBE2110-20-1.5	Comb, 20 sample, 1.5 mm thick
SBE2110-CP	Mini cooling pack


Electroblotting Accessories (for use with 10 x 10 system)

Catalog No.	Description
SBE2010-BM	Blotting module, includes 3 cassettes
SBE2110-B-MC	Mini blotting cassette
SBE2110-B-FP	Fiber pads for blotting cassette, pk of 6
SBE2110-CP	Mini cooling pack

2D Electrophoresis Accessories (tube gel accessories only)

Catalog No.	Description
SBE2010-2M	Tube gel module with blanking ports
SBE2110-2D-CT	Capillary tubes, disposable, pk of 100
SBE2110-2D-BP	Blanking ports, pk of 10
SBE2110-CP	Mini cooling pack

PAGE
Module



Tube Gel
Module



Electroblotting
Module



Select
BioProducts
For Life Science Research

Email: sales@selectbioproducts.com
www.selectbioproducts.com