

Force Mini



The Force Mini is a personal centrifuge that is ideal for quick spin downs of microtubes and PCR tubes. The centrifuge was designed to be extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm's. Press the lid release button to open the lid and the rotor quickly decelerates to stop for removing samples.

- Ideal for microfiltration and quick spins
- Complete with 8 x 2.0mL microtube rotor and 4 x 8 PCR strip rotor
- Quickly reaches maximum speed of 6,000 rpm/2,000 x g
- Transparent lid for easy sample viewing

Ordering Information:

Cat. No.	Description
SBC140	Force Mini with 2.0mL tube and strip tube rotors, eight adapters for 0.4 and 0.5mL tubes, 115V
SBC140-2	Force Mini with 1.5mL tube and strip tube rotors, eight adapters for 0.4 and 0.5mL tubes, 230V EU
SBC140-3	Force Mini with 1.5mL tube and strip tube rotors, eight adapters for 0.4 and 0.5mL tubes, 230V UK
SBC140TR	Additional 8 x 1.5mL tube rotor for Force Mini
SBC140SR	Additional strip tube rotor for Force Mini
C1205	Adapters for 0.5/0.6mL tubes, pack of 6
C1206	Adapters for 0.25/0.4mL tubes, pack of 6
C1222	Adapters for 0.2mL tubes, pack of 6

